

ARIES[®]

ARIES[®] Systems | Welcome to the New Way to Work

ARIES[®] *Bordetella* Assay



The ARIES[®] *Bordetella* Assay is a real-time polymerase chain reaction (PCR) based qualitative in vitro diagnostic test for the direct detection and identification of *Bordetella pertussis* and *Bordetella parapertussis* nucleic acid in nasopharyngeal swab (NPS) specimens obtained from individuals suspected of having a respiratory tract infection attributable to *B. pertussis* or *B. parapertussis*. Features include:

- **Specific:** Aids in diagnosis and surveillance reporting by identifying patients with *B. pertussis* infection
- **Flexible and Efficient:** Run 1 to 12 tests per batch, utilizing both STAT testing and low to medium sample batching
- **Fully Integrated:** Automate all aspects of testing, from sample preparation through analysis
- **Fast Time to Results:** Answers in less than 2 hours with minimal hands-on time results
- **Error-Reducing Safeguards:** Internal barcode scanning matches samples to cassettes and may reduce data input errors

Performance

The performance of the moderate complexity ARIES[®] *Bordetella* Assay was assessed at five (5) geographically diverse clinical sites in the United States.

Refer to Package Insert for additional details: Luminex Corporation | ARIES[®] *Bordetella* Assay (IVD) Kit Package Insert.

Table 1: ARIES[®] *Bordetella* Assay Performance for *B. pertussis*

Specimen Description	PPA		95% CI	NPA		95% CI
Prospective	30/32*	93.8%	79.2% - 99.2%	1009/1020	98.9%	98.1% - 99.5%
Pre-selected	37/37	100%	90.5% - 100%	77/77	100%	95.3% - 100%
Total	67/69	97.1%	89.9% - 99.6%	1086/1097	99.0%	98.2% - 99.5%


*Two (2) prospective specimens generated false negative results by ARIES[®] *Bordetella* Assay when compared to the composite comparator method (02-179 and 06-267).

Table 2: ARIES[®] *Bordetella* Assay Performance for *B. parapertussis*

Specimen Description	PPA		95% CI	NPA		95% CI
Prospective	2/2	100%	15.8% - 100%	1048/1050	99.8%	99.3% - 100%
Pre-selected	20/20	100%	83.2% - 100%	93/94*	98.9%	94.2% - 100%
Contrived	50/50	100%	92.9% - 100%	50/50	100%	92.9% - 100%
Total	72/72	100%	95.0% - 100%	1191/1194	99.7%	99.3% - 99.9%

*One (1) pre-selected specimen generated a false positive result by ARIES[®] *Bordetella* Assay when compared to the composite comparator method (01-122).

Workflow



The operator simply adds specimen to the sample chamber, puts the cassette in the magazine, loads the magazine into the ARIES® System, and the run will start automatically.

Usage and Targets

The direct detection and identification of *B. pertussis* and *B. parapertussis* nucleic acids from symptomatic patients aids in the diagnosis of *B. pertussis* and *B. parapertussis* respiratory infection in conjunction with other clinical findings and epidemiological information. Negative results for the ARIES® *Bordetella* Assay do not preclude *B. pertussis* or *B. parapertussis* infection and positive results do not rule out co-infections with other respiratory pathogens. The ARIES® *Bordetella* Assay targets the *B. pertussis* toxin promoter and the *B. parapertussis* IS1001 insertion element in the bacterial genomes. The ARIES® *Bordetella* Assay is indicated for use with ARIES® Systems.

Ordering Information*

Product Name	Part Number
ARIES® <i>Bordetella</i> Assay	50-10037 (24 tests)
ARIES® <i>Bordetella</i> Assay Protocol File	CN-0388-01 (one time order only)
ARIES® Two Module System Includes: Instrument System Operation Manual Two Magazines	Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand ARIES-M12V1-IVD
ARIES® M1 System Includes: Instrument System Operation Manual One Magazine	Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand ARIES-M6V1-IVD
SYNCT™ Software	CN-SW47

*Products are CE Marked and FDA Cleared for IVD use.

Luminex®
complexity simplified.

orders@luminexcorp.com or support@luminexcorp.com

For In Vitro Diagnostic Use. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. Validation of the LIS compatibility must be performed by the end user. ARIES® Systems are class 1(I) laser products.

©2017 Luminex Corporation. All rights reserved. Luminex and ARIES are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

www.luminexcorp.com/bordetella

HEADQUARTERS

UNITED STATES
+1.512.219.8020
info@luminexcorp.com

EUROPE
+31.73.800.1900
europe@luminexcorp.com

CANADA
+1.416.593.4323
info@luminexcorp.com

CHINA
+86.21.802.31150
info@luminexcorp.com

JAPAN
+81.3.5545.7440
info@luminexcorp.com

SS26437