



NxTAG[®] Respiratory Pathogen Panel v2 (CA-IVD)

The NxTAG[®] Respiratory Pathogen Panel (RPP) v2 is a multiplex, real-time polymerase chain reaction-based test, that qualitatively detects nucleic acid from 2 bacterial and 20 viral pathogens, including severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), from nasopharyngeal swab (NPS) samples in universal transport medium (UTM[™]), MicroTest[™] M4RT[™], or viral transport media (VTM).

Panel-based respiratory testing results can help to reduce hospitalization time and unnecessary or prolonged antibiotic courses.¹ While respiratory illnesses often present with similar symptoms, panel-based assays like NxTAG RPP v2 offer an advantage to clinicians looking to identify the causative pathogen(s) underlying a patient's symptoms. With 22 respiratory targets differentiated and scalable throughput, included and scalable throughput, you can count on this assay to support testing even when sample volumes surge.

Designed for use on the NxTAG[®]-Enabled MAGPIX[®] System, this assay was designed to keep pace with virus evolution and is the enhanced, second-generation version (v2) of the original NxTAG RPP.

1. Making Sense of Respiratory Viral Panel Results. ASM (Internet). Cited February 2023. Available from: <https://asm.org/Articles/2020/March/Making-Sense-of-Respiratory-Viral-Panel-Results>.

The NxTAG[®] Respiratory Pathogen Panel v2 Assay offers

- **Comprehensive Detection of Respiratory Pathogens:** Detects 2 bacterial and 20 viral pathogens, including SARS-CoV-2 (ORF1ab and M gene) in a single test, enabling accurate diagnosis and patient treatment.
- **Scalable Throughput:** Process up to 96 samples in less than 3 hours post-extraction, accommodating variable, day-to-day testing demands.
- **Minimal Hands-On Time:** Pre-plated, lyophilized reagents facilitate a simple workflow with just one pipetting step, ensuring an easy fit in any lab's daily routine.

Performance

Limit of Detection (LoD) for SARS-CoV-2

The development of version 2 focused on overall improved performance and the addition of the SARS-CoV-2 targets. The LoD for SARS-CoV-2 is 500 copies/mL (Strain: USA-WA 1/2020).

Targets

| Viral Targets | | Bacterial Targets |
|-------------------------------|----------------------------|----------------------------------|
| Adenovirus | Influenza A subtype H3 | <i>Chlamydomphila pneumoniae</i> |
| Coronavirus 229E | Influenza B | <i>Mycoplasma pneumoniae</i> |
| Coronavirus HKU1 | Parainfluenza 1 | |
| Coronavirus NL63 | Parainfluenza 2 | |
| Coronavirus OC43 | Parainfluenza 3 | |
| Human Bocavirus | Parainfluenza 4 | |
| Human Metapneumovirus | Rhinovirus/ Enterovirus | |
| Influenza A | RSV A | |
| Influenza A subtype H1 | RSV B | |
| Influenza A subtype 2009 H1N1 | SARS-CoV-2 | |

Summary of the Comparative Clinical Performance of the NxTAG® RPP v2 Assay

| Pathogen | PPA* | | NPA* | |
|----------------------------------|--------------|-------|--------------|-------|
| | TP / (TP+FN) | % | TN / (TN+FP) | % |
| Virus | | | | |
| Adenovirus | 31/33 | 93.9% | 1261/1264 | 99.8% |
| Coronavirus 229E | 12/12 | 100% | 1284/1285 | 99.9% |
| Coronavirus HKU1 | 30/32 | 93.8% | 1265/1265 | 100% |
| Coronavirus NL63 | 30/30 | 100% | 1268/1268 | 100% |
| Coronavirus OC43 | 30/30 | 100% | 1263/1267 | 99.7% |
| Human bocavirus | 8/8 | 100% | 1416/1419 | 99.8% |
| Human Metapneumovirus | 30/30 | 90% | 1265/1267 | 99.8% |
| Influenza A | 39/39 | 100% | 1257/1258 | 99.9% |
| Influenza A subtype H1 | 0 | N/A | 1296/1297 | 99.9% |
| Influenza A subtype 2009 H1N1 | 29/30 | 96.7% | 1267/1267 | 100% |
| Influenza A subtype H3 | 9/9 | 100% | 1288/1288 | 100% |
| Influenza B | 30/30 | 100% | 1267/1267 | 100% |
| Parainfluenza 1 | 30/30 | 100% | 1266/1267 | 99.9% |
| Parainfluenza 2 | 30/30 | 100% | 1267/1267 | 100% |
| Parainfluenza 3 | 35/35 | 100% | 1262/1262 | 100% |
| Parainfluenza 4 | 15/16 | 93.8% | 1281/1281 | 100% |
| Rhinovirus/Enterovirus | 112/115 | 97.4% | 1155/1182 | 97.7% |
| RSV A | 30/30 | 100% | 1265/1267 | 99.8% |
| RSV B | 0 | N/A | 1297/1297 | 100% |
| SARS-CoV-2 | 65/66 | 98.5% | 1354/1361 | 99.5% |
| Bacteria | | | | |
| <i>Chlamydomphila pneumoniae</i> | 14/14 | 100% | 1282/1283 | 99.9% |
| <i>Mycoplasma pneumoniae</i> | 48/52 | 92.3% | 1245/1245 | 100% |

*Additional details are available in the Package Insert.

Influenza A and B inclusivity was assessed with sequences available from the GISAID database between January 1, 2017 and January 6, 2022. The assay oligos for influenza A, influenza A H1, influenza A H3, and influenza B are predicted to have 100% inclusivity against the analyzed sequences.

We are continuously monitoring the GISAID database SARS-CoV-2 Variants of Concern sequences. Please go to <https://www.luminexcorp.com/covid19-testing-solutions/> and select "Click here" to view the results of our *in silico* inclusivity analysis.

Workflow



Ordering Information

| Product Name | Part Number |
|--|----------------------|
| NxTAG® Respiratory Pathogen Panel v2 CA-IVD (96 tests) | I055C0465 |
| NxTAG®-Enabled MAGPIX® System | MAGPIX-XPON4.1-CEIVD |

Luminex®
A DiaSorin Company

orders@luminexcorp.com or support@luminexcorp.com

©2023 Luminex Corporation. A DiaSorin Company. All rights reserved. NXTAG and MAGPIX are trademarks of Luminex Corporation, registered in the US and other countries. All other trademarks are property of their respective owners.

[luminexcorp.com/NxTAG-RPP/](https://www.luminexcorp.com/NxTAG-RPP/)

HEADQUARTERS

UNITED STATES
+1 512 219 8020
info@luminexcorp.com

EUROPE
+31 73 800 1900
europe@luminexcorp.com

CANADA
+1 416 593 4323
info@luminexcorp.com

CHINA
+86 21 8036 9888
info@luminexcorp.com

JAPAN
+81 3 5545 7440
info@luminexcorp.com

SS535205