



ARIES[®]

Flu A/B & RSV Assay (RUO)



Welcome to the New Way to Work


The ARIES[®] Flu A/B & RSV Assay (RUO)* is a rapid, accurate method for the detection and differentiation of three key respiratory pathogens. Research laboratories can identify influenza A virus, influenza B virus, and respiratory syncytial virus (RSV) from nasopharyngeal swab (NPS) specimens using a sample to answer platform. Key features include:

- Integrated sample processing control ensures the assay run is successful from extraction through amplification
- Internal barcode scanning matches samples to cassettes, enabling Position Independent Results no matter where each cassette is placed
- Bidirectional LIS connectivity enables electronic ordering and reduces data input errors
- Simultaneous STAT and batch testing on a single ARIES[®] System enables reduced turnaround times
- Ability to mask targets empowers the user to report only the results of interest



The ARIES[®] System is crafted to increase laboratory efficiency, ensure result accuracy, and fit seamlessly into today's lean laboratory.

Workflow



The operator simply adds specimen to the sample chamber, puts the cassette in the magazine, loads the magazine into the ARIES® System, and the run will start automatically.

Ordering Information

Description	Part Number
ARIES® Flu A/B & RSV Assay Kit (RUO)*	50-10014 (24 tests)
ARIES® Flu A/B & RSV Assay Protocol File (RUO)*	CN-0327-01 (one time order only)
ARIES® Two Module System (RUO)* Includes: Instrument System Operation Manual Two Magazines Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner	ARIES-M12V1-RUO

Luminex®

support@luminexcorp.com or orders@luminexcorp.com

* RUO - For Research Use Only. Not for use in diagnostic procedures. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. The ARIES® System is a class 1(I) laser product. Validation of the LIS compatibility must be performed by the end user.

©2016 Luminex Corporation. All rights reserved. Luminex and ARIES are trademarks of Luminex Corporation, registered in the U.S. and other countries.

HEADQUARTERS

UNITED STATES
+1.512.219.8020
info@luminexcorp.com
www.luminexcorp.com

EUROPE
+31.73.800.1900
europe@luminexcorp.com
www.luminexcorp.com

CANADA
+1.416.593.4323
info@luminexcorp.com
www.luminexcorp.com

CHINA
+86.21.802.31150
chinainfo@luminexcorp.com
www.luminexcorp.com

JAPAN
+81.3.5545.7440
infojp@luminexcorp.com
www.luminexcorp.com

551485.0216