



Certificate of Quality

MAGPIX® Calibration Kit, IVD MPXIVD-CAL-K25

| | | |
|---|---------------------|----------------|
| REF | Product Number | MPXIVD-CAL-K25 |
| LOT | Kit Lot Number | |
|  | Kit Expiration Date | |
|  | Date of Manufacture | |

| Kit Components | Part Number | Lot* |
|--|-----------------|------|
| MAGPIX Calibrator, Microspheres, 6 mL | MPXCAL-05 | |
| 8-Well Microtiter plate strips, Pack of 25 | 13-52047 | N/A |
| MAGPIX Calibration Kit CD (contains Kit Lot Importable File, component CoQs, and Kit Instructions) | 89-20286-00-001 | |

This certificate provides descriptive lot information and assurance that quality requirements have been met for the components, intermediates, and finished product. Luminex Corporation certifies that the MAGPIX Calibration Kit, IVD, MPXIVD-CAL-K25 lot defined above conforms to all prescribed acceptance criteria and is fit for its intended use.

*Reagent target values may vary lot to lot. Only the reagent lots listed should be used when performing system calibration with this kit lot.

Luminex®

EN ISO 13485 Certified Company

12212 Technology Blvd • Austin, Texas 78727-6115 USA

Phone: 1.512.381.4397 • Toll-free: 1.877.785.2323 • Fax: 512.219.5195 • www.luminexcorp.com