

ARIES[®]

ARIES[®] Systems

ARIES[®] Group A Strep Assay (RUO)



The ARIES[®] Group A Strep Assay is a rapid method for the direct detection of *Streptococcus pyogenes* (Group A β -hemolytic *Streptococcus*) nucleic acid from throat swab samples using a sample to answer platform.

- **Flexibility** — Run as many as 12 tests per batch, supporting STAT testing and low- to medium-volume batch testing
- **Fast Time to Results** — Get answers in less than two hours with minimal hands-on time
- **Full Integration** — Automate all aspects of testing, from sample extraction through analysis, to enable easy implementation across multiple shifts
- **Error-reducing Safeguards** — Internal barcode scanning matches samples to cassettes and may reduce data input errors



ARIES[®] increases research laboratory efficiency, ensures result accuracy, and fits seamlessly into today's lean laboratory.

Workflow



The operator prepares the sample, adds it to the cassette, puts the cassette in the magazine, loads the magazine into the ARIES® System, and the run will initiate on its own if the auto-run feature has been enabled.

Ordering Information

Product Name	Part Number
ARIES® Group A Strep Assay (RUO)	50-10040 (24 tests)
ARIES® Group A Strep Assay Protocol File Kit (RUO)	CN-0390-01 (one time order only)
ARIES® Two Module System Includes: Instrument System Operation Manual Two Magazines	Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand ARIES-M12V1-RUO
ARIES® M1 System Includes: Instrument System Operation Manual One Magazine	Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand ARIES-M6V1-RUO
SYNCT™ Software	CN-SW34

Luminex®
complexity simplified.

orders@luminexcorp.com or support@luminexcorp.com

This product is for Research Use Only. Not for use in diagnostic procedures. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. ARIES® Systems are class 1(I) laser products. Validation of the LIS compatibility must be performed by the end user.

©2017-2019 Luminex Corporation. All rights reserved. Luminex and ARIES are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

www.luminexcorp.com/gas-ruo/

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europe@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@cn.luminexcorp.com

JAPAN

+81.3.5545.7440

info@jp.luminexcorp.com

SS121920