

Luminex[®]

ARIES[®]

ARIES[®] Systems | Welcome to the New Way to Work

ARIES[®] *Bordetella* Assay (RUO)




The ARIES[®] *Bordetella* Assay (RUO) is a rapid, accurate method for the detection and differentiation of two key *Bordetella* species. With this assay, research laboratories can identify *B. pertussis* and *B. parapertussis* from specimens such as nasopharyngeal swabs (NPS), using a sample to answer platform. Key features include:

- **Full Integration** — Automate all aspects of testing, from sample extraction through analysis, to enable easy implementation across multiple shifts
- **Flexibility** — Run as many as 12 tests per batch, supporting STAT testing and low- to medium-volume batch testing
- **Error-reducing Safeguards** — Internal barcode scanning matches samples to cassettes and may reduce data input errors



ARIES[®] Systems are crafted to increase laboratory efficiency, ensure result accuracy, and fit seamlessly into today's lean laboratory.

Workflow



The operator simply adds specimen to the sample chamber, puts the cassette in the magazine, loads the magazine into an ARIES® System, and the run will start automatically.

Ordering Information

Description	Part Number
ARIES® Bordetella Assay Kit	50-10036 (24 tests)
ARIES® Bordetella Assay Protocol File	CN-0387-01 (one time order only)
ARIES® Two Module System Includes: Instrument System Operation Manual Two Magazines Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand	ARIES-M12V1-RUO
ARIES® M1 System Includes: Instrument System Operation Manual One Magazine Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand	ARIES-M6V1-RUO
SYNCT™ Software	CN-SW34

Luminex®
complexity simplified.

orders@luminexcorp.com or support@luminexcorp.com

For Research Use Only. Not for use in diagnostic procedures. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. The ARIES® System is a class 1(I) laser product. Validation of LIS compatibility must be performed by the end user.

©2018 Luminex Corporation. All rights reserved. Luminex and ARIES are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

www.luminexcorp.com/bordetella-ruo/

HEADQUARTERS
UNITED STATES
+1.512.219.8020
info@luminexcorp.com

EUROPE
+31.73.800.1900
europe@luminexcorp.com

CANADA
+1.416.593.4323
info@luminexcorp.com

CHINA
+86.21.8036.9888
info.cn@luminexcorp.com

JAPAN
+81.3.5545.7440
info.jp@luminexcorp.com

SS121953