

**Luminex**<sup>®</sup>

# ARIES<sup>®</sup>

ARIES<sup>®</sup> Systems

## ARIES<sup>®</sup> GBS Assay



The ARIES<sup>®</sup> GBS Assay is a rapid method for the detection of *Streptococcus agalactiae*, also known as Group B *Streptococcus*, from LIM broth enriched specimens using a sample to answer platform. The ARIES<sup>®</sup> GBS Assay (RUO) offers:

- **Fast Time to Results** — Get answers in less than two hours with minimal hands-on time
- **Flexibility** — Run as many as 12 tests per batch, supporting STAT testing and low- to medium-volume batch testing
- **Full Integration** — Automate all aspects of testing, from sample extraction through analysis, to enable easy implementation across multiple shifts
- **Error-reducing Safeguards** — Internal barcode scanning matches samples to cassettes and may reduce data input errors



ARIES<sup>®</sup> Systems are crafted to increase research laboratory efficiency, ensure result accuracy, and fit seamlessly into today's lean laboratory.

## Workflow



The operator simply adds specimen to the sample chamber, puts the cassette in the magazine, loads the magazine into the ARIES® System, and the run will initiate on its own if the auto-run feature has been enabled.

## Ordering Information

Product Name	Part Number
ARIES® GBS Assay (RUO)	50-10015 (24 tests)
ARIES® GBS Assay Protocol File (RUO)	CN-0328-01 (one time order only)
ARIES® Two Module System Includes: Instrument System Operation Manual Two Magazines	Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand ARIES-M12V1-RUO
ARIES® M1 System Includes: Instrument System Operation Manual One Magazine	Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand ARIES-M6V1-RUO
SYNCT™ Software	CN-SW34

**Luminex®**  
complexity simplified.

**orders@luminexcorp.com or support@luminexcorp.com**

This product is for Research Use Only. Not for use in diagnostic procedures. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. ARIES® Systems are class 1(l) laser products. Validation of the LIS compatibility must be performed by the end user.

©2015-2018 Luminex Corporation. All rights reserved. Luminex and ARIES are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

[www.luminexcorp.com/gbs-ruo/](http://www.luminexcorp.com/gbs-ruo/)

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europe@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@luminexcorp.com

JAPAN

+81.3.5545.7440

infojp@luminexcorp.com

SS122291