



NxTAG[®]

Respiratory Pathogen Panel



Syndromic Testing Made Easy


The NxTAG[®] Respiratory Pathogen Panel detects 21 respiratory pathogens in an easy to use, single well, closed tube system.

- **Excellent Clinical Performance:** Gain confidence in your results with excellent clinical performance
- **Comprehensive Panel:** Easily detect 21 important pathogens in a single test
- **Rapid Turnaround Time:** Decrease your hands-on time and turnaround time with a quick, easy protocol
- **Scalable Throughput:** Scalable throughput lets you run 1-96 samples at a time
- **Customizable:** Select any number of targets within the panel to customize the assay

Viral Targets		
Influenza A	Parainfluenza 1	Human Rhinovirus / Enterovirus
Influenza A H1	Parainfluenza 2	Adenovirus
Influenza A H3	Parainfluenza 3	Coronavirus HKU1
Influenza B	Parainfluenza 4	Coronavirus NL63
Respiratory Syncytial Virus A (RSV A)	Human Bocavirus	Coronavirus OC43
Respiratory Syncytial Virus B (RSV B)	Human Metapneumovirus (hMPV)	Coronavirus 229E
Bacterial Targets		
<i>Chlamydomphila pneumoniae</i>	<i>Legionella pneumophila</i>	<i>Mycoplasma pneumoniae</i>

Save Time and Streamline Your Workflow

Get Reliable Data with Minimal Hands-on Time



Time estimated on 24 samples run simultaneously.

Simple and Speedy Workflow

Step 1

Add 1-96 extracted samples to pre-plated test wells



Step 2

Integrated multiplex PCR and bead hybridization



Step 3

Read on MAGPIX[®]





To learn more, consult with your local Luminex representative.

<https://www.luminexcorp.com/NxTAG-RPP>

For In Vitro Diagnostic Use. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence.

©2015 Luminex Corporation. All rights reserved. The trademarks mentioned herein are the property of Luminex or their respective owners.

HEADQUARTERS

UNITED STATES
 +1.512.219.8020
info@luminexcorp.com
www.luminexcorp.com

EUROPE
 +31.73.800.1900
europa@luminexcorp.com
www.luminexcorp.eu

CANADA
 +1.416.593.4323
info@luminexcorp.com
www.luminexcorp.com

CHINA
 +86.21.802.31150
chinainfo@luminexcorp.com
www.luminexcorp.com

JAPAN
 +81.3.5545.7440
infojp@luminexcorp.com
www.luminexcorp.com