

NxTAG®

For use with the MAGPIX® System

NxTAG® Respiratory Pathogen Panel



Broad Panel Testing Made Easy

The NxTAG® Respiratory Pathogen Panel detects 20 respiratory pathogens in an easy to use, single well, closed tube system.

The panel's unique features and benefits include:

- **Excellent Performance:** Have confidence in the results and reduce validation burden¹
- **Comprehensive Content:** Detect 20 important pathogens in a single test
- **Rapid Turnaround Time:** Decrease your hands-on time and turnaround time with a quick, easy protocol
- **Scalable Throughput:** Run from 1-96 samples at a time
- **Customizable Panel:** Select any number of targets within the assay

Viral Targets


- Influenza A
- Influenza A H1
- Influenza A H3
- Influenza B
- Respiratory Syncytial Virus A
- Respiratory Syncytial Virus B
- Parainfluenza virus 1
- Parainfluenza virus 2
- Parainfluenza virus 3
- Parainfluenza virus 4
- Human Bocavirus
- Human Metapneumovirus
- Rhinovirus/Enterovirus
- Adenovirus
- Coronavirus HKU1
- Coronavirus NL63
- Coronavirus OC43
- Coronavirus 229E

Bacterial Targets

- Chlamydomphila pneumoniae*
- Mycoplasma pneumoniae*

Save Time and Streamline Your Workflow


Get Reliable Data with Minimal Hands-on Time



■ NxTAG Time to Results ■ NxTAG Hands-on Time

Time estimated on 24 samples run simultaneously

Workflow



Step 1

Add 1-96 extracted samples to pre-plated test wells



Step 2

Integrated multiplex PCR and bead hybridization



Step 3

Read on MAGPIX®

Ordering Information

Product Name	Part Number
NxTAG® Respiratory Pathogen Panel (US-IVD)	I051C0447 (96 tests)
MAGPIX® System	MAGPIX-XPON4.1-CE-IVD

The results of this test should not be used as the sole basis for diagnosis, treatment, or other patient management decisions. Please refer to the IVD package insert for the full intended use, limitations, and risk information.

REFERENCES

1. Luminex Corporation | NxTAG Respiratory Pathogen Panel US-IVD Package Insert.

Luminex
complexity simplified.

orders@luminexcorp.com or support@luminexcorp.com

For In Vitro Diagnostic Use. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence.

©2016-2019 Luminex Corporation. All rights reserved. Luminex, NxTAG, and MAGPIX are trademarks of Luminex Corporation, registered in the U.S. and other countries.

www.luminexcorp.com/NxTAG

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europa@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@luminexcorp.com

JAPAN

+81.3.5545.7440

info@luminexcorp.com

SS147782