

SYNCT™

SYNCT™ Software | The Power Behind Your Lab


SYNCT™ Standard Curve Analysis Software SYNCT QUANT 2.0



SYNCT™ Standard Curve Analysis Software (SYNCT Quant 2.0) is specifically created to enable standard curve analysis for Laboratory Developed Tests (LDTs). Labs can create quantitative or qualitative LDTs using the powerful ARIES® UDP App and with this latest version of SYNCT Quant 2.0, standard curve creation occurs within the application. LDTs can also have a standard curve assigned to them so quantitative results are automatically calculated at the end of the run without any necessary user interaction. Features and benefits include:

- Ability to create standard curves within the application, eliminating the need to use external tools for curve creation
- Automatically apply standard curve to results for faster reporting turnaround
- Bidirectional LIS interface, enabling the transfer of both orders and results, both qualitative and quantitative
- One connection between SYNCT Software and the LIS, allowing labs to easily connect multiple instruments
- Results monitoring and management from anywhere in the lab for faster reporting times
- Ability to associate cassettes to orders from anywhere to allow for workflow customization and optimization

Using SYNCT Quant 2.0 to create quantitative results from your LDTs can be done using these simple steps:



Ordering Information

Part Number	Item
CN-SW59-02	SYNCT™ Standard Curve Analysis Software (SYNCT™ Quant 2.0)
CN-SW39-01	ARIES® UDP App*

*SYNCT is required to run the ARIES® UDP (User Defined Protocol) App.

WARNING: This version of SYNCT is For Laboratory Use and cannot process IVD or RUO assay results. Use of IVD or RUO assays may require separate appropriately-labeled versions of SYNCT and ARIES® Systems.

Luminex®
complexity simplified.

orders@luminexcorp.com or support@luminexcorp.com

For Laboratory Use. For use with SYNCT Software and ARIES® Systems. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. ARIES® Systems are class 1(I) laser products.

©2019 Luminex Corporation. All rights reserved. ARIES and Luminex are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

luminexcorp.com/SYNCT/

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europe@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@cn.luminexcorp.com

JAPAN

+81.3.5545.7440

info@jp.luminexcorp.com

SS163142