


SYNCT[™]

SYNCT[™] Software | The Power Behind Your Lab

SYNCT[™] Standard Curve Analysis Software




SYNCT[™] Standard Curve Analysis Software (SYNCT 1.1u3) is software specifically created to enable standard curve analysis for Laboratory Developed Tests (LDTs). Labs can still create LDTs using the powerful ARIES[®] UDP App, but now you will have the added functionality of saving and applying standard curves to calculate quantitative results.

In addition to the capability to manage orders from Laboratory Information Systems (LIS), view test results, prepare test reports, perform administrative functions, and maintain security control, this software offers:

- The ability to save standard curves to easily create quantitative LDT results
- A bidirectional LIS interface, enabling the transfer of orders and results, both qualitative and quantitative
- One connection between SYNCT Software and the LIS, allowing labs to easily connect multiple instruments
- Results monitoring and management from anywhere in the lab
- The ability to associate cassettes to orders from anywhere where testing is done

Using SYNCT 1.1u3 to create quantitative results from your LDTs can be done with these simple steps:



Ordering Information

Part Number	Item
CN-SW59-01	SYNCT™ Standard Curve Analysis Software (SYNCT 1.1u3)*
CN-SW52-01	ARIES® UDP App*

*SYNCT is required to run the ARIES® UDP (User Defined Protocol) App.

WARNING: This version of SYNCT is For Laboratory Use and cannot process IVD or RUO assay results. Use of IVD or RUO assays may require separate appropriately-labeled versions of SYNCT and ARIES® Systems.

Luminex
complexity simplified.

orders@luminexcorp.com or support@luminexcorp.com

For Laboratory Use. For use with SYNCT Software and ARIES® Systems. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence. ARIES® Systems are class 1(I) laser products.

©2018 Luminex Corporation. All rights reserved. ARIES and Luminex are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

www.luminexcorp.com/SYNCT/

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europe@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@luminexcorp.com

JAPAN

+81.3.5545.7440

infojp@luminexcorp.com

SS93707